

# URGOPORE



## PRODUCT PRESENTATION

### On plastic spool without cover

CODE	DIMENSIONS	MINIMUM DELIVERY UNIT
1795	9.14 m x 1.25 cm	80 rolls
1796	9.14 m x 2 cm	60 rolls
1797	9.14 m x 2.5 cm	50 rolls
1798	9.14 m x 5 cm	30 rolls
1799	9.14 m x 7.5 cm	10 rolls
1785	5 m x 1.25 cm	80 rolls
1786	5 m x 2 cm	60 rolls
1787	5 m x 2.5 cm	50 rolls
1788	5 m x 5 cm	30 rolls
1789	5 m x 7.5 cm	10 rolls

### On plastic spool with cover

CODE	DIMENSIONS	MINIMUM DELIVERY UNIT
1772	5 m x 1.25 cm	30 rolls
1774	5 m x 2 cm	30 rolls
1775	5 m x 2.5 cm	30 rolls
1776	5 m x 5 cm	10 rolls
1781	5 m x 7.5 cm	10 rolls
1915	8 m x 2.5 cm dispenser	200 rolls

NON-WOVEN  
MICROPOROUS  
ADHESIVE TAPE  
HIGH CUTANEOUS  
TOLERANCE

# URGOPORE

## PRODUCT DESCRIPTION

- Supple non-woven microporous backing.
- High skin tolerance acrylic adhesive composition.
- X-ray transparent.

## PRODUCT ADVANTAGES

- Micro-ventilated : permits skin breathing and elimination of exudation and transpiration liquids from the wound thus avoiding maceration.
- High skin tolerance : minimizes risks on the most sensitive skin.
- Painless removal.

## INDICATIONS

- All departments involving fixing and maintaining of dressings and medical material requiring a frequent renewal on allergic, fragile or devitalized skin.

