

URGOSYVAL



PRODUCT PRESENTATION

On plastic spool without cover

CODE	DIMENSIONS	MINIMUM DELIVERY UNIT
1765	5 m x 1.25 cm	80 rolls
1746	5 m x 2 cm	60 rolls
1747	5 m x 2.5 cm	50 rolls
1748	5 m x 5 cm	30 rolls
1769	5 m x 7.5 cm	30 rolls

On plastic spool with cover

CODE	DIMENSIONS	MINIMUM DELIVERY UNIT
1741	5 m x 1.25 cm	30 rolls
1742	5 m x 2 cm	30 rolls
1743	5 m x 2.5 cm	30 rolls
1744	5 m x 5 cm	10 rolls
1757	5 m x 7.5 cm	10 rolls

Perforated in individual box

CODE	DIMENSIONS	MINIMUM DELIVERY UNIT
1758	5 m x 10 cm	10 rolls
1759	5 m x 18 cm	10 rolls

HIGH SKIN
TOLERANCE
TAFFETA ADHESIVE
TAPE

URGOSYVAL

PRODUCT DESCRIPTION

- White taffeta textile backing.
- High skin tolerance solvent-free adhesive composition.

PRODUCT ADVANTAGES

- Cellulose acetate backing which can be cut in both directions without scissors.
- High adhesive power, good skin tolerance due to its composition of purified components.
- Large sizes aerated through perforations.

INDICATIONS

- For non-occlusive but long lasting dressings or for skin sensitive to terpenic resins or colophon.
- Application of medical device on sensitive skins.

