


# URGOSTRIPS



**PRODUCT PRESENTATION**

**Box of 50 foils :**

CODE	DIMENSIONS	NUMBER OF SKIN CLOSURES PER FOIL
1600	75 mm x 3 mm	5
1601	75 mm x 6 mm	3
1602	100 mm x 6 mm	10
1603	100 mm x 12.5 mm	6

## STERILE ADHESIVE SKIN CLOSURES

# URGOSTRIPS

## PRODUCT DESCRIPTION

- Non-woven polyester backing, aerated, set on a silicon paper. Tear resistance =  $4.9 \pm 1$  kg/cm.
- Acrylic adhesive mass with high cutaneous tolerance.
- Sterilised with ethylene oxide.
- Wrapped in an individual peel-off sachet transparent on one side.

## PRODUCT ADVANTAGES

- The aeration of the microporous backing avoid maceration.
- Security of use :
  - . High tear resistance obtained without adding reinforcement threads.
  - . High and long-lasting adhesiveness without risk of sliding and respect of the different skin types.
- Ease of use :
  - . Double self-parting strip which enables prehension of the product in respect with aseptic action.
  - . Peel-off sachet.
- Easy, accelerated and aesthetic healing.
- Painless application and removal.

## INDICATIONS

- Repairer surgery :
  - . mastectomy,
  - . mending of scar in post-surgery.
- Dermatological surgery : cutaneous tumour's ablation (Nevus...)
- Emergencies : quick and efficient joining of the edges of the wound.
- To reinforce the traditional stitches.

